

ACCESSORIES

Swing out rotor 6 place, treated steel.



Part No. **7001391**

Buckets and adapters with engraved identification. Autoclavable.



Buckets with two compartments for blood bags of 500 to 800 ml. For triple and quad bags
Part No. **7001392**

Single bucket. Capacity: 1 blood bag of 500 to 800 ml or 1 flask of 1 litre or with tube adapters.
Part No. **7001393**



Adapters for single bucket Part No. 7001393.

Bottle



1001394



7001395



7001396



7001397



7001398



7001399



7001400



7001401

Double Bucket	Single Bucket	Part No. adapter	Total max. Volume (6 adapters)	Capacity: bags, tubes by adapter	Ø x Length mm	R.p.m. Min./Max.	R.c.f. (xg) Min./Max.
7001392	-	-	9600 ml	2 bags 500 or 800 ml	-	300/4070	29/5370
-		Bottle 1000 ml 1001394	6000 ml	1 bottle 1000 ml or 1 bag	98 x 175	300/4070	29/5370
-	7001393	7001395	3000 ml	1 bottle 500 ml	70 x 165	300/4070	29/5370
-	Capacity	7001396	1200 ml	1 bottle 200 ml	56.6 x 104	300/4070	29/5370
-	without adapter:	7001397	1200 ml	2 tubes 100 ml	44 x 104	300/4070	29/5370
-	1 bag of	7001398	1200 ml	4 tubes 50 ml	34 x 100	300/4070	29/5370
-	500 to 800 ml	7001399	1350 ml	9 tubes 25 ml	24 x 100	300/4070	29/5370
-	or 1x bottle	7001400	1530 ml	17 tubes 15 ml	16.5 x 100	300/4070	29/5370
-	1000 ml.	7001401	1008 ml	24 tubes 7 ml	12 x 100	300/4070	29/5370

Note: All buckets and adapters come without tubes unless stated, see page 127 for tubes and accessories.

Other adapters can be manufactured to the customers specification.