


24




Coloana cromatograf cu sicane si slif mama robinet cu dop de sticla

Chromatography column with indentations socket and stopcock with glass plug

COD CODE	CAPACITATE CAPACITY	DIAMETRU INTERIOR INSIDE DIAMETER	LUNGIME UMLERE FILL LENGTH
	mL	mm	mm
23.160.02	15	10	200
23.160.04	120	20	400
23.160.06	430	30	600
23.160.08	1000	40	800

25



Coloana cromatograf cu sicane si slif mama robinet cu dop de teflon

Chromatography column with indentations socket and stopcock with PTFE plug

COD CODE	CAPACITATE CAPACITY	DIAMETRU INTERIOR INSIDE DIAMETER	LUNGIME UMLERE FILL LENGTH
	mL	mm	mm
23.162.02	15	10	200
23.162.04	120	20	400
23.162.06	430	30	600
23.162.08	1000	40	800

26




Coloana cromatograf cu placa filtranta de porozitate 0, robinet cu dop de sticla

Chromatography column with fused in filter disc porosity 0 and stopcock with glass plug

COD CODE	CAPACITATE CAPACITY	DIAMETRU INTERIOR INSIDE DIAMETER	LUNGIME UMLERE FILL LENGTH
	mL	mm	mm
23.164.02	15	10	200
23.164.04	35	15	200
23.164.06	70	15	400
23.164.08	125	20	400
23.164.10	190	20	600
23.164.12	280	30	400
23.164.14	430	30	600
23.164.16	1000	40	800

27



Coloana cromatograf cu placa filtranta de porozitate 0, robinet cu dop de teflon

Chromatography column with fused in filter disc porosity 0 and stopcock with PTFE plug

COD CODE	CAPACITATE CAPACITY	DIAMETRU INTERIOR INSIDE DIAMETER	LUNGIME UMLERE FILL LENGTH
	mL	mm	mm
23.166.02	15	10	200
23.166.04	35	15	200
23.166.06	70	15	400
23.166.08	125	20	400
23.166.10	190	20	600
23.166.12	280	30	400
23.166.14	430	30	600
23.166.16	1000	40	800




34


Coloana Hempel fara manta

Hempel column with packing support, with cone and socket without jacket

COD CODE	LUNGIME EFECTIVA EFFECTIVE LENGTH	DIAMETRU INTERIOR INSIDE DIAMETER	SLIF TATA / MAMA CONE NS / SOCKET NS
	mm	mm	
24.172.02	200	15	14/23
24.172.04	300	15	14/23
24.172.06	300	30	29/32
24.172.08	500	30	29/32


35


Coloana Vigreux fara manta

Vigreux column with cone and socket, indentations alternately horizontal and sloping, without jacket

COD CODE	LUNGIME EFECTIVA EFFECTIVE LENGTH	DIAMETRU INTERIOR INSIDE DIAMETER	SLIF TATA / MAMA CONE NS / SOCKET NS
	mm	mm	
24.128.02	200	25	14/23
24.128.04	300	25	14/23
24.128.01	200	25	29/32
24.128.03	350	25	29/32
24.128.08	200	25	29/32
24.128.10	300	25	29/32
24.128.11	500	25	29/32
24.128.12	600	25	29/32


23



Coloana cromatograf cu sicane robinet cu dop de teflon
Chromatography column with indentations and stopcock with PTFE plug

COD CODE	CAPACITATE CAPACITY	DIAMETRU INTERIOR INSIDE DIAMETER	LUNGIME UMPLERE FILL LENGTH
	mL	mm	mm
23.158.02	15	10	200
23.158.04	35	15	200
23.158.06	70	15	400
23.158.08	125	20	400
23.158.10	190	20	600
23.158.12	280	30	400
23.158.14	430	30	600
23.158.16	1000	40	800

22



Coloana cromatograf cu sicane robinet cu dop de sticla
Chromatography column with indentations and stopcock with glass plug

COD CODE	CAPACITATE CAPACITY	DIAMETRU INTERIOR INSIDE DIAMETER	LUNGIME UMPLERE FILL LENGTH
	mL	mm	mm
23.156.02	15	10	200
23.156.04	35	15	200
23.156.06	70	15	400
23.156.08	125	20	400
23.156.10	190	20	600
23.156.12	280	30	400
23.156.14	430	30	600
23.156.16	1000	40	800