

Eye protection / safety goggles



Ultra light build of only 18 grammes total weight. Maximum flexibility, perfect and slip-free fit.

Safety goggles Ultralight

Color	Price €	Quantity scale from	Scale price €	Item no.
blue		5		2502-1002 N3



Ultra-light protective goggles with sporty design, maximum wearing comfort, according to DIN/EN standard.

Casual style safety goggles

Color	Price €	Quantity scale from	Scale price €	Item no.
blue/orange		5		2502-2010 N3



Over-goggles with hingeless side pieces and with soft side piece ends. Combines the advantages of goggles with side pieces with the wearing properties of headbands. The flexible side piece design allows it to adapt to different wearing situations, with or without glasses. DIN EN 166 - 168 and 170.

Safety goggles Ultraflex

Color	Price €	Quantity scale from	Scale price €	Item no.
blue		5		2502-1003 N3



With anatomically designed wrap-around screen for an increased degree of safety. Unrestricted view in all directions, larger cover by 20% compared with normal glasses with side pieces, UV protection up to 400 nm.

Safety goggles Style

Color	Price €	Quantity scale from	Scale price €	Item no.
black		5		2502-1006 N3



Lightweight, unrestricted lateral sight, variable side pieces and side piece length which can be set to 4 different lengths guarantee the perfect fit without pressure points due to Duo-Flex side pieces, lateral safety shields integrated in the goggle's lens, according to DIN EN 166-168 and 170.

Astrolite safety goggles

Color	Price €	Quantity scale from	Scale price €	Item no.
black/blue		5		2502-4000 N3
black/red		5		2502-4001 N3

Eye protection / safety goggles



Anti-fog lens, suitable for wearers of glasses, with ventilation system, unrestricted lateral view, close fitting, pressure-free fit, according to DIN/EN standard.



Various colours, pleasantly light. Good adjustability from the option for setting the side pieces to 4 lengths. Shield made of polycarbonate, resistant to high impact with integrated lateral protection. Maximum UV protection.



Very comfortable over-goggles with unrestricted side view, suitable to go over most glasses. 180° field of vision, lateral ventilation through the side shields. Side pieces can be adjusted to 4 lengths. Without coating. With Polycarbonate side pieces. According to DIN EN 166 -168 and 170.

► Ideal as goggles for visitors



Extremely robust eye-wear case with reinforced zipper and belt loop. Suitable for all spectacles, not for safety goggles Ultravision.



Very robust large eye-wear case for panoramic goggles with zipper, belt loop and additional compartment. Suitable for all spectacles and panoramic goggles.

Ultravision safety goggles

Color	Price €	Quantity scale from	Scale price €	Item no.
orange		5		2502-1007 N3

Color safety spectacles

Color	Price €	Quantity scale from	Scale price €	Item no.
yellow/black		5		2502-1004 N3
blue		5		2502-1005 N3
silver gray		5		2502-1010 N3
violet/green		5		2502-1012 N3

Panorama safety goggles

Color	Price €	Quantity scale from	Scale price €	Item no.
glass clear		5		2502-1001 N3

Case with belt loop for all glasses with temple

Color	Price €	Quantity scale from	Scale price €	Item no.
black		5		2502-1008 N3

Case with belt loop for all glasses and goggles

Color	Price €	Quantity scale from	Scale price €	Item no.
black		5		2502-1009 N3

Plastic apron made of PVC

With eyelets, neck band and waist band welded into the apron, industry quality, black, fabric reinforced, size 80 cm x 100 cm.

! Attention! Make sure to pay attention to the "List of Chemical Resistance"!



Plastic apron made of PVC

Price €	Quantity scale from	Scale price €	Item no.
	10		2507-1001 N3

SafetyBox



! Attention! Make sure to pay attention to the "List of Chemical Resistance"!

Panorama safety goggles and protective gloves for personal safety, stored cleanly in a box. Recommended by trade associations.

SafetyBox

Price €	Quantity scale from	Scale price €	Item no.
	10		2507-9000 N3

Drain safety stand



Residual liquid from pump tubes or samplers can be reliably collected with the drain safety stand. The safety stand is resistant to diluted bases and acids. The collecting container can be filled with a neutralisation medium which dilutes and neutralises residual liquid.

The stand pipe is mounted on the wall with sturdy pipe clamps. The pipe clamps can simply be opened to remove the safety stand.

The container for collecting liquid can be unscrewed, making it convenient to replace or empty.

Drain safety stand

Price €	Quantity scale from	Scale price €	Item no.
	5		5618-1000 N4