


19



Refrigerent Liebig cu slif tata si mama DIN 12567
Liebig-condenser with cone and socket DIN 12567

COD CODE	SLIF TATA / MAMA	LUNGIME MANTA REFRIGERENT
	CONE NS/ SOCKET NS	JACKET LENGTH
		mm
24.140.02	14/23	160
24.140.04	14/23	200
24.140.06	14/23	300
24.140.08	29/32	200
24.140.10	29/32	250
24.140.12	29/32	300
24.140.14	29/32	400
24.140.16	29/32	500


20



Refrigerent Allihn DIN 12581 cu slif mama si tata
Allihn-condenser with cone and socket DIN 12581

COD CODE	SLIF TATA / MAMA	LUNGIME MANTA REFRIGERENT
	CONE NS/ SOCKET NS	JACKET LENGTH
		mm
24.142.02	14/23	160
24.142.04	14/23	200
24.142.06	14/23	300
24.142.08	29/32	200
24.142.10	29/32	250
24.142.12	29/32	300
24.142.14	29/32	400
24.142.16	29/32	600
24.142.07	19/26	300
24.142.15	29/32	220
24.142.17	45/40	300


21



Refrigerent cu serpentina,cu slif mama si tata
Spiral condenser with cone and socket

COD CODE	SLIF TATA / MAMA	LUNGIME MANTA REFRIGERENT
	CONE NS/ SOCKET NS	JACKET LENGTH
		mm
24.144.02	14/23	160
24.144.04	14/23	200
24.144.06	29/32	250
24.144.08	29/32	300
24.144.10	29/32	400


22



Refrigerent Dimroth cu slif tata
Dimroth-condenser with cone

COD CODE	SLIF TATA / MAMA	LUNGIME MANTA REFRIGERENT
	CONE NS/ SOCKET NS	JACKET LENGTH
		mm
24.146.02	14/23	200
24.146.04	29/32	200
24.146.06	29/32	300
24.146.08	34/35	250
24.146.10	45/40	300


23



Refrigerent Dimroth cu slif mama si slif tata DIN 12591
Dimroth-condenser with double coil and cone DIN 12591

COD CODE	SLIF TATA / MAMA	LUNGIME MANTA REFRIGERENT
	CONE NS/ SOCKET NS	JACKET LENGTH
		mm
24.150.02	14/23	160
24.150.04	14/23	200
24.150.06	14/23	300
24.150.08	29/32	200
24.150.10	29/32	250
24.150.12	29/32	300
24.150.14	29/32	400


24



Refrigerent Dimroth cu manta dubla cu slif mama si tata
Jacketed coil condenser with cooling jacket and cooling coil

COD CODE	SLIF TATA / MAMA	LUNGIME MANTA REFRIGERENT
	CONE NS/ SOCKET NS	JACKET LENGTH
		mm
24.152.02	14/23	200
24.152.04	29/32	200
24.152.06	29/32	250
24.152.08	29/32	300
24.152.10	29/32	400

18



Refrigerent Liebig cu slif tata
Liebig-condenser with cone

COD CODE	SLIF TATA	LUNGIMEA TUBULUI
	CONE NS	TUBE LENGTH
	mm	mm
24.138.02	14/23	160
24.138.04	14/23	200
24.138.06	14/23	300
24.138.08	29/32	200
24.138.10	29/32	250
24.138.12	29/32	300
24.138.14	29/32	400


17



Refrigerent de aer cu slif tata
Air condenser with cone


COD CODE	SLIF TATA	LUNGIMEA TUBULUI
	CONE NS	TUBE LENGTH
		mm
24.136.02	14/23	700
24.136.04	*29/32	1000



*
1
4

Refrigerent Dimroth cu spirala dubla
Dimroth-condenser with double spiral

COD CODE	SLIF TATA CONE NS	LUNGIME MANTA REFRIGERENT TUBE LENGTH	
25.108.02	45/40	300	mm
25.108.04	60/46	300	

*
1
3

Refrigerent Dimroth
Dimroth-condenser

COD CODE	SLIF TATA CONE NS	LUNGIME MANTA REFRIGERENT TUBE LENGTH	
			mm
24.146.04	29/32	200	
24.146.08	34/35	250	
24.146.10	45/40	300	