

































## Angle Rotors (High Speed Centrifuge)

	Tube Capacity	0.2ml	0.5ml	1.5ml	2.0ml	<b>Mega 17R</b> Max. RPM : 15,000 Max. RCF : 29,104 xg <b>Supra 22K</b> Max. RPM : 18,000 Max. RCF : 41,910 xg <b>Supra 30K</b> Max. RPM : 20,000 Max. RCF : 51,741 xg
	<b>A1.5S-36</b> No.14 DIM. Ø x L mm	6 x 18	8 x 30	11 x 38	11 x 38	
Tube Type						
Adapter						
Tube per Rotor	36					
Type	Angle Rotor ∠45°C	Max.RCF	42,573	46,643	52,546	
Max. RPM	20,000 rpm	Radius(mm)	95.2	104.3	115.7	
Max. RCF	51,741xg	Boring Ø x L mm	6 x 21	8 x 37	11 x 37.5	
Radius	115.7 mm					
	Tube Capacity	10ml				<b>Mega 17R</b> Max. RPM : 17,000 Max. RCF : 27,496 xg <b>Supra 22K</b> Max. RPM : 22,000 Max. RCF : 46,049 xg <b>Supra 30K</b> Max. RPM : 30,000 Max. RCF : 85,628 xg
	<b>A10S-12</b> No.2 DIM. Ø x L mm	16 x 82				
Tube Type						
Adapter						
Tube per Rotor	12					
Type	Angle Rotor ∠36°C	Max.RCF	85,628			
Max. RPM	30,000 rpm	Radius(mm)	85.1			
Max. RCF	85,628 xg	Boring Ø x L mm	16.3 x 74.5			
Radius	85.1 mm					
	Tube Capacity	15ml				<b>Mega 17R</b> Max. RPM : 15,000 Max. RCF : 25,004 xg <b>Supra 22K</b> Max. RPM : 17,000 Max. RCF : 32,116 xg <b>Supra 30K</b> Max. RPM : 17,000 Max. RCF : 32,116 xg
	<b>A15S-12C</b> No.8 DIM. Ø x L mm	17 x 120				
Tube Type						
Adapter						
Tube per Rotor	12					
Type	Angle Rotor ∠25°C	Max.RCF	32,116			
Max. RPM	17,000 rpm	Radius(mm)	99.4			
Max. RCF	32,116 x g	Boring Ø x L mm	17 x 115			
Radius	99.4 mm					
	Tube Capacity	15(16)ml				<b>Mega 17R</b> Max. RPM : 17,000 Max. RCF : 31,050 xg <b>Supra 22K</b> Max. RPM : 22,000 Max. RCF : 52,001 xg <b>Supra 30K</b> Max. RPM : 30,000 Max. RCF : 96,696 xg
	<b>A15S-12</b> No.4 DIM. Ø x L mm	18 x 100				
Tube Type						
Adapter						
Tube per Rotor	12					
Type	Angle Rotor ∠25°C	Max.RCF	96,696			
Max. RPM	30,000 rpm	Radius(mm)	96.1			
Max. RCF	96,696 xg	Boring Ø x L mm	19 x 96			
Radius	96.1 mm					
	Tube Capacity	50ml				<b>Mega 17R</b> Max. RPM : 15,000 Max. RCF : 22,891 xg <b>Supra 22K</b> Max. RPM : 17,000 Max. RCF : 29,402 xg <b>Supra 30K</b> Max. RPM : 17,000 Max. RCF : 29,402 xg
	<b>A50S-6C</b> No.6 DIM. Ø x L mm	29 x 115				
Tube Type						
Adapter						
Tube per Rotor	6					
Type	Angle Rotor ∠40°C	Max.RCF	29,402			
Max. RPM	17,000 rpm	Radius(mm)	91.0			
Max. RCF	29,402 xg	Boring Ø x L mm	29.5 x 108.6			
Radius	91 mm					

## Angle Rotors (High Speed Centrifuge)

<b>A50S-6 No.5</b> 	Tu Tube Capacity	15(16)ml	50ml	50ml	<b>Mega 17R</b> Max. RPM : 17,000 Max. RCF : 31,050 xg <b>Supra 22K</b> Max. RPM : 22,000 Max. RCF : 52,001 <b>Supra 30K</b> Max. RPM : 30,000 Max. RCF : 96,696 xg
	DIM. Ø x L mm	18 x 100	28.7 x 104	28.8 x 107	
Tube Type					
Adapter					
Tube per Rotor	6				
Max.RCF	90,558	96,696			
ra <sub>1</sub> Radius(mm)	90.0	96			
bo Boring x L mm	18.5 x 92	29 x 100			
Ca.....					
Type Angle Rotor ∠30°C					
Max. RPM 30,000 rpm					
Max. RC <sup>RCF</sup> 96,696 xg					
Radius 96.1 mm					
<b>A50S-8 No.7</b> 	Tu Tube Capacity	15(16)ml	50ml	50ml	<b>Mega 17R</b> Max. RPM : 15,000 Max. RCF : 24,878 xg <b>Supra 22K</b> Max. RPM : 22,000 Max. RCF : 53,516 xg <b>Supra 30K</b> Max. RPM : 30,000 Max. RCF : 99,513 xg
	DIM. Ø x L mm	18 x 100	28.7 x 104	28.8 x 107	
Tube Type					
Adapter					
Tube per Rotor	8				
Max.RCF	93,375	99,513			
ra <sub>1</sub> Radius(mm)	92.8	98.9			
bo Boring x L mm	18.5 x 92	29 x 100			
Ca.....					
Type Angle Rotor ∠36°C					
Max. RPM 30,000 rpm					
Max. RC <sup>RCF</sup> 99,513 xg					
Radius 98.9 mm					
<b>A85S-6 No.9</b> 	Tu Tube Capacity	85ml			<b>Mega 17R</b> Max. RPM : 15,000 Max. RCF : 23,193 xg <b>Supra 22K</b> Max. RPM : 20,000 Max. RCF : 41,232 xg <b>Supra 30K</b> Max. RPM : 25,000 Max. RCF : 64,425 xg
	DIM. Ø x L mm	38.2 x 105.7			
Tube Type					
Adapter					
Tube per Rotor	6				
Max.RCF	64,425				
ra <sub>1</sub> Radius(mm)	92.2				
bo Boring x L mm	38.5 x 90				
Ca.....					
Type Angle Rotor ∠25°C					
Max. RPM 25,000 rpm					
Max. RC <sup>RCF</sup> 64,425 xg					
Radius 92.2 mm					
<b>A250S-6 No.10</b> 	Tu Tube Capacity	250ml	250ml		<b>Supra 22K</b> Max. RPM : 12,000 Max. RCF : 22,265 xg <b>Supra 30K</b> Max. RPM : 14,000 Max. RCF : 30,305 xg
	DIM. Ø x L mm	61.8 x 127.7	61 x 134		
Tube Type					
Adapter					
Tube per Rotor	6				
Max.RCF	30,305				
ra <sub>1</sub> Radius(mm)	138.3				
bo Boring x L mm	62 x 103				
Ca.....					
Type Angle Rotor ∠25°C					
Max. RPM 14,000 rpm					
Max. RC <sup>RCF</sup> 30,305 xg					
Radius 138.3mm					