

Antibiotic test paper ref. 1468

Recommended for resistance test of pathogen organisms



APLICATIONS

- Impregnation of antibiotics for resistance test of pathogen organisms.

Absorbent paper circles without impregnation and very pure. Provided in 6, 9 and 13mm diameter for its use in antibiotic test.

Formats and measures



6	9	13
---	---	----

Dimensions: Diameter en mm
Presentation: Boxes of 1.000 units



20 x 20	40 x 40	60 x 60
---------	---------	---------

Dimensions: Measure in cm
Presentation: Packs of 100 units

Other formats and dimensions available under demand.

Information for orders. Antibiotic test paper ref. 1468

Format	Diameter (mm)	Cod.	units/pack
Circles	6	PN1468006	1.000
Circles	9	PN1468009	1.000
Circles	13	PN1468013	1.000

Format	Size: width x length (cm)	Cod.	units/pack
Sheets	20 x 20	HJ14682020	100
Sheets	40 x 40	HJ14684040	100
Sheets	60 x 60	HJ14686060	100

Multi Lab® SRL Romania
www.multilab.ro
contact@multilab.ro service@multilab.ro info@multilab.ro
Rel.Client 0730 711101 **Hot-Line 0724044800**
Desfacere 072 8930531 / 073 5876690 / 0735876689 / 0728930530
Service 0728930535 / 0728930536 / 0731309088
Secretariat (+04) 021 3127144 / 45 Fax. (+04) 021 3127141