

OVERHEAD STIRRERS

Accesorii pentru agitatoare mecanice FALC



Stainless steel stirring shafts

Arbori de amestecare din otel inoxidabil

Code	Photo n°	Rotor mm length x H mm	Ø rod mm	Overall length mm
605.0152.02	1	50x10	7	300
605.0152.03	2	50x20	7	300
605.0152.04	3	70x50	7	300
605.0152.05	4	45x40	8	350
605.0152.06	5	100x10	7	300
605.0152.07	1	100x20	7	300
605.0152.08	6	50	7	300
605.0152.10	4	90x80	8	350

SS stirring shafts with blade PTFE

Arbori de amestecare din otel inoxidabil cu lama din PTFE

Code	Rotor mm length x H mm	Ø rod mm	PTFE thickness	Overall length mm
605.0152.12	90x20	6,5	1,5	450
605.0152.13	90x20	6,5	1,5	600



SS stirring shafts with one paddle PTFE

Arbor de amestecare cu o singura paleta din PTFE

Code	Rotor mm length x H mm	Ø rod mm	PTFE thickness	Overall length mm
605.0152.15	80x12	6,5	1,5	450
605.0152.16	80x12	6,5	1,5	600



SS U-shaped stirrer shafts PTFE

Arbori de amestecare cu paleta forma U din PTFE

Code	Rotor mm length x H mm	Ø rod mm	PTFE thickness	Overall length mm
605.0152.17	60x40	6,5	1,5	350
605.0152.18	60x40	6,5	1,5	450
605.0152.19	80x50	6,5	1,5	450



SS fan shaped stirrer shafts PTFE

Arbori de amestecare cu paleta forma Vane din PTFE

Code	Rotor mm length x H mm	Ø rod mm	PTFE thickness	Overall length mm
605.0152.23	38	6,5	1,5	300
605.0152.22	24	6,5	1,5	300



Supports e clamps Stative si cleme de prindere

Code	Description
263.8350.20	Support opening 130x280mm with one rod ØxH12x600
263.8350.22	Support opening 180x460mm with one rod ØxH15x750
263.8350.72	Support opening 230x600mm with three rods ØxH 15x1000mm
146.2900.41	Claw clamp with maximum Ø18mm for rods

