


SWING-OUT ROTOR RS 4/750

SMART LOCK - rotor exchange without any tool
Up to 4 x 750 ml centrifuge tubes



Rotor specifications	RS 4/750
Max. speed	4.500 rpm
Max. G force	4.530 rcf
Radius	200 mm
Angle	90°
Max. capacity	4 x 750 ml
Min. acceleration time	38 s
Min. braking time	31 s
Rotor autoclavable	yes

Code	234.000
Order ref.	424161



Buckets (set of 4)	B 750
Code	234.030
Order ref.	424162
Buckets autoclavable	yes



Lids for B 750(set of 4)	BL 750
Code	234.050
Order ref.	409269
Lid autoclavable	no



MTP Carriers (set of 4) - for 4 standard or 1 deep well plate	BM 750
Code	233.600
Order ref.	424151
Buckets autoclavable	yes

	Blood collection tubes	Blood collection tubes BD	Blood collection tubes Terumo	Blood collection tubes	Conical tubes	Conical tubes
Nominal volume	1,0 - 5 ml	1,0 - 5 ml	1,0 - 5 ml	7 - 12 ml	15 ml	50 ml
Tube dimensions (d x l)	13 x 75 mm 13 x 100 mm	13 x 75 mm 13 x 100 mm	13 x 75 mm 13 x 100 mm	16 x 100 mm	17 x 120 mm	30 x 115 mm
Inserts (set of 4)	I 750 22/5-7	I 750 27/5-7	I 750 31/5-7	I 750 20/10	I 750 15/15	I 750 7/50
No. of tubes / insert	22	27	31	20	15	7
No. of tubes / rotor	88	108	124	80	60	28
Radius	196 mm	198 mm	196 mm	198 mm	198 mm	168 mm
G Force	4.440 rcf	4.480 rcf	4.440 rcf	4.480 rcf	4.480 rcf	3.800 rcf
Code	234.060	234.080	234.100	234.120	234.140	234.160
Order ref.	424163	424164	424165	424166	424167	424168
Inserts autoclavable	yes	yes	yes	yes	yes	yes