

ACCESSORIES

Swing out rotor with interchangeable sealed buckets and adapters, each with identification code. Autoclavable.



Part No. head	Rotor type	Bucket open	Sealed bucket	Part No. adapter	Total max. Volume (4 adapters)	Bucket/Adapter capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg)
	Swing out	7000731	7000732	Bottle 200 ml 5715202**	800 ml	1 x 200 ml	56.6 x 104	5000	4137
	Swing out	7000731	7000732	7000733	400 ml	1 x 100 ml	44 x 104	5000	4137
	Swing out	7000731	7000732	7000734	200 ml	1 x 50 ml	34 x 100	5000	4137
	Swing out	7000731	7000732	7000935	200 ml	1 x 50 ml Falcon	30 x 117	5000	4137
	Swing out	7000731	7000732	7001645	240 ml	2 x 30 ml	25 x 92	5000	4137
	Swing out	7000731	7000732	7000735	200 ml	2 x 25 ml	24 x 100	5000	4137
	Swing out	7000731	7000732	7000736	420 ml	7 x 15 ml and VAC*	16.5 x 100	5000	4137
7009079	Swing out	7000731	7000732	7000934	180 ml	3 x 15 ml Falcon	16.5 x 120	5000	4137
	Swing out	7000731	7000732	7001646	-	4 x 17 Ø x 110 mm	-	5000	4137
	Swing out	7000731	7000732	7001647	240 ml	6 x 10 ml	17 x 110	5000	4137
	Swing out	7000731	7000732	7000737	336 ml	12 x 7 ml	12 x 100	5000	4137
	Swing out	7000731	7000732	7000738	240 ml	12 x 5 ml	12.5 x 75	5000	4137
	Swing out	7000731	7000732	7001648	140 ml	7 x 5 ml and VAC*	13 x 75	5000	4137
	Swing out	7000731	7000732	7000739	144 ml	12 x 3 ml	11 x 100	5000	4137
	Swing out	7000731	7000732	7000740	-	undrilled	-	5000	4137
7000609	Swing out	-	-	-	-	6 x 96 microtiter	-	4200	2702

*VAC: Blood tube with stopper.

** Bucket part nos. 7000731 and 7000732 accommodate 1 x 200 ml polypropylene bottle with cap without the need of an adapter.

Note: All buckets and adapters come without tubes unless stated, see page 132 for tubes and accessories.

HEAD REDUCERS

Reducers for head. Part No. 7000741:

- 7001770 For 1 tube 85 ml (38 Ø x 105 mm).
- 7001771 For 1 tube 50 ml (34 Ø x 100 mm).
- 7001772 For 1 tube 50 ml (29 Ø x 103 mm).
- 7001773 For 1 tube 25 ml (24 Ø x 100 mm).
- 7001774 For 1 tube 30 ml (25 Ø x 100 mm).
- 7001775 For 3 tubes 15 ml (16,5 Ø x 100 mm).
- 7001776 For 3 tubes 10 ml (16 Ø x 80 mm).
- 7001777 For 7 tubes 7 ml (12 Ø x 100 mm).
- 7001778 For 7 tubes 5 ml (12,5 Ø x 75 mm).

Reducers for head. Part No. 7000742:

- 7001707 For 1 tube 50 ml (29 Ø x 103 mm).
- 7001708 For 1 tube 25 ml (24 Ø x 100 mm).
- 7001709 For 1 tube 30 ml (25 Ø x 92 mm).
- 7001710 For 1 tube 15 ml (16,5 Ø x 100 mm).

7001780 For 1 tube 10 ml (16 Ø x 80 mm).

- 7001781 For 4 tubes 7 ml (12 Ø x 100 mm).
- 7001782 For 4 tubes 5 ml (12,5 Ø x 75 mm).

Reducers for head. Part No. 7001449:

- 7001783 For 1 tube 50 ml (34 Ø x 100 mm).
- 7001784 For 1 tube 50 ml (29 Ø x 103 mm).
- 7001785 For 1 tube 25 ml (24 Ø x 100 mm).
- 7001786 For 1 tube 30 ml (25 Ø x 92 mm).
- 7001787 For 2 tubes 15 ml (16,5 Ø x 100 mm).
- 7001788 For 2 tubes 10 ml (16 Ø x 100 mm).
- 7001789 For 5 tubes 7 ml (12 Ø x 100 mm).
- 7001790 For 5 tubes 5 ml (12,5 Ø x 75 mm).

Reducers for head. Part No. 7001448:

- 7001791 For 1 tube 25 ml (24 Ø x 100 mm).
- 7001792 For 1 tube 30 ml (25 Ø x 100 mm).

7001793 For 1 tube 15 ml (16,5 Ø x 100 mm).

- 7001794 For 1 tube 10 ml (16 Ø x 80 mm).
- 7001795 For 2 tubes 7 ml (12 Ø x 100 mm).
- 7001796 For 2 tubes 5 ml (12,5 Ø x 75 mm).

Reducers for head. Part No. 7001447:

- 7001797 For 1 tube 15 ml (16,5 Ø x 100 mm).
- 7001798 For 1 tube 10 ml (16 Ø x 80 mm).
- 7001799 For 1 tube 7 ml (12 Ø x 100 mm).
- 7001800 For 1 tube 5 ml (12,5 Ø x 75 mm).

Reducers for head. Part No. 7001445/7001446:

- 7001801 For 1 tube 5 ml (12,5 Ø x 75 mm).

Reducers for head. Part No. 7000743:

- 7000700 For 1 tube 7 ml (12 Ø x 100 mm).
- 1001557 For 1 tube 5 ml (12,5 Ø x 75 mm).