

COMMON ACCESSORIES. Tubes of glass and plastic.

Plastic tubes and flasks with screw cap lids.

Volume ml	Dimensions mm	Type	Part No.	Material	Max. (xg) range	Temperature °C
10	16 Ø x 80	Comecta	5752399	Polypropylene	50000	4 ... 121
		Herolab	5752421	Polypropylene	57000	-180 ... 145
		Herolab	5752401	Polycarbonate	57000	-150 ... 121
30	25 Ø x 92	Comecta	5752398	Polypropylene	50000	4 ... 121
		Herolab	5752955	Polypropylene	64000	-180 ... 145
		Herolab	5752951	Polycarbonate	64000	-150 ... 121
50	29 Ø x 103	Comecta	5752397	Polypropylene	50000	4 ... 121
		Herolab	5753205	Polypropylene	69000	-180 ... 145
		Herolab	5753201	Polycarbonate	69000	-150 ... 121
80	38 Ø x 105	Comecta	5752396	Polypropylene	40000	4 ... 121
		Herolab	5753243	Polypropylene	50000	-180 ... 145
		Herolab	5753241	Polycarbonate	50000	-150 ... 121
200	56.6 Ø x 104	Herolab	5715202	Polypropylene	10000	-180 ... 145
		Herolab	5715203	Polycarbonate	10000	-150 ... 121
250	62 Ø x 122	Comecta	5752395	Polypropylene	13200	4 ... 121
		Herolab	5753411	Polypropylene	33000	-180 ... 145
		Herolab	5753401	Polycarbonate	33000	-150 ... 121
500	70 Ø x 165	Herolab	5753512	Polypropylene	14000	-180 ... 145
		Herolab	5753501	Polycarbonate	14000	-150 ... 121
1000	98 Ø x 175	Herolab	5752571	Polypropylene	4936	-180 ... 145
		Herolab	5752561	Polycarbonate	4936	-150 ... 121



Type HEROLAB polycarbonate tubes transparent Polypropylene tubes translucent. Autoclavable.

Type COMECTA polypropylene translucent. Autoclavable.



Glass Tubes.

Volume	Dimensions mm	Part No.
Micro-haematocrit*	1.4 Ø x 75	7000910
7 ml	12 Ø x 100	7000773
15 ml	16.5 Ø x 100	7000772
25 ml	24 Ø x 100	7000771
50 ml	34 Ø x 100	7000770
100 ml	44 Ø x 104	7000769
250 ml	56 Ø x 147	7000768
ASTM	Conical	7001088
ASTM	Pear shape	7000767

*Box of 100 tubes.

SPARE

Buckets for angle rotors.

Suitable for model	Part No. 7001403 15 ml plastic	Part No. 7001500 15 ml plastic	Part No. 7005000 50 ml metal	Part No. 7001000 100 ml metal
Cencom II Part No. 7002240	Yes			
Centro-8 Part No. 7002356		Yes		
Centromix BLT Part No. 7002556		Rotor 7000550		
Meditronic BLT Part No. 7001476		Rotor 7000743	Rotor 7000742	Rotor 7000741
Medifriger BLT Part No. 7001477		Rotor 7000743	Rotor 7000742	Rotor 7000741
Macrotronic BLT Part No. 7001485		Rotor 7000560	Rotor 7000561	
Macrofriger BLT Part No. 7001486		Rotor 7000560	Rotor 7000561	

ALL THE ROTORS, BUCKETS AND ADAPTERS COME WITH LASER ETCHED IDENTIFICATION.

ACCESSORIES for cytology

Description	Part No.
1. Slide clip fastener	5004050
2. Sample container	5004051
3. Two samples container	5004052
4. Filter	5004053
5. Cyto slide of 1 circle	5004054
6. Cyto slide of 2 circles	5004055
7. Disposable container 1 sample	5004056
8. Disposable container 2 samples	5004057

