

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Angle rotor	12×15ml/10/5ml	6000r/min	6080xg	PP/PC tube
	Adaptor	12×2ml			PP/PC tube
NO.2 	Swing rotor	steel frame rotor	5000r/min	4752xg	PP/PC tube
	Swing frame	4×100ml			PP/PC tube
	Adaptor	4×15/10ml			PP/PC tube
	Swing frame	4×50ml			PP/PC tube
	Adaptor	4×15/10ml			PP/PC tube
NO.3 	Swing rotor	steel rotor	4000r/min	3041xg	PP/PC tube
	Swing frame	4×2×100ml			PP/PC tube
	Adaptor	4×2×15/10ml			PP/PC tube
	Swing frame	4×2×50ml			PP/PC tube
	Adaptor	4×2×15/10ml			PP/PC tube
	Swing frame	4×8×15/10/5ml			PP/PC tube
	Rack	4×12×5/2ml			Vacuum blood collection tubes
NO.4 	Swing rotor	4X18×5/2ml	4000r/min	2773xg	Vacuum blood collection tubes
NO.5 	Swing bucket	4×600ml	4000r/min	3354xg	Square cup with wind cover
	Adapter	4×25×5/2ml			Vacuum blood collection tubes
	Adapter	4×12×15/10ml			PP/PC tube
	Adapter	4×5×50ml			PP/PC tube
	Adapter	4×4×100ml			PP/PC tube
NO.6 	Swing bucket	6×600ml	4000r/min	3400xg	Square cup with wind cover
	Adapter	6×25×5/2ml			Vacuum blood collection tubes
	Adapter	6×12×15/10ml			PP/PC tube
	Adapter	6×5×50ml			PP/PCtube
	Adapter	6×4×100ml			PP/PCtube

NO.7 	Swing bucket	4×500ml	4000r/min	3354xg	Round cup with wind cover
	Adapter	4×17×5/2ml			Vacuum blood collection tubes
	Adapter	4×8×15/10ml			PP/PC tube
	Adapter	4×3×50ml			PP/PC tube
	Adapter	4×100ml			PP/PC tube
NO.8 	Swing bucket	6×500ml	4500r/min	4324xg	Round cup with wind cover
	Adapter	6×17×5/2ml			Vacuum blood collection tubes
	Adaptor	6×8×15/10ml			PP/PC tube
	Adapter	6×3×50ml			PP/PC tube
	Adapter	6×100ml			PP/PC tube
NO.9 	Swing bucket	4×800ml	4000r/min	3443xg	Square cup with wind cover
	Adapter	4×30×5/2ml			Vacuum blood collection tubes
	Adapter	4×16×15/10ml			PP/PC tube
	Adapter	4×6×50ml			PP/PC tube
	Adapter	4×4×100ml			PP/PC tube
NO.10 	Swing bucket	4×750ml	4500r/min	4437xg	Round cup
	Adapter	4×26×5/2ml			Vacuum blood collection tubes
	Adapter	4×13×15/10ml			PP/PC tube
	Adapter	4×5×50ml			PP/PC tube
	Adapter	4×3×100ml			PP/PC tube
NO.11 	Swing bucket	4×1000ml	4000r/min	3578xg	Rhombus cup with wind cover
	Adapter	4×42×5/2ml			Vacuum blood collection tubes
	Adapter	4×23×15/10ml			PP/PC tube
	Adapter	4×7×50ml			PP/PC tube
	Adapter	4×4×100ml			PP/PC tube
NO.12 	Microplate rotor	4×2×96hole	4200r/min	2859xg	Square cup 4*2*96 deep well plate 4*4*96 shallow well plate
	Adapter	4×32×1.5ml		3057xg	PP/PC tube