

**APPLICATION**

Wall mounting

**DETAILS**

Waterflow around 13l/min  
at 3 Bars pressure

Test pressure: 9 Bars/air

Max working pressure: 10 Bars

**CONTACT MATERIALS**

PP

Ceramic

**HEADWORK SEALING**

CERAMIC temperature  
working range: 0-50°

Headwork in PP with CERAMIC  
sealing without oil, grease  
and lubricants

**CONNECTION TAP**

Male thread G1/2"

J.G. for tube Ø10

**STANDARD AND NORMS**

EN13792:2000

DIN12898


DIN12918

ISO228/1

**OPTIONS**

Straight J.G. connector  
11010.2MDSPP

90° J.G. connector  
11010.2MDSPPB



**APPLICATION**

Wall mounting

**DETAILS**

Waterflow around 13l/min  
at 3 Bars pressure

Test pressure: 9 Bars/air

Max working pressure: 10 Bars

**CONTACT MATERIALS**

PP

Ceramic

**HEADWORK SEALING**

CERAMIC temperature  
working range: 0-50°

Headwork in PP with CERAMIC  
sealing without oil, grease  
and lubricants

**CONNECTION TAP**

Male thread G1/2"

J.G. for tube Ø10

**STANDARD AND NORMS**

EN13792:2000

DIN12898

DIN12918

ISO228/1

**OPTIONS**

Straight J.G. connector  
11012.2MDSPP

90° J.G. connector  
11012.2MDSPPB

