

STEREO Series - Accessories

ACCESSORIES FOR STX AND MS-2

ST-009	Dust cover.
ST-015	Glass stage, 60 mm diameter.
M-114	C-Mount adapter for 1/2" sensor.
M-115	C-Mount adapter for 1/3" sensor.
M-116	C-Mount adapter for 2/3" sensor.
M-173	Photo adapter for APS-C and Full Frame Reflex cameras.
M-005	Micrometer slide for software calibration, 1mm/10um, 10mm/100um.
15104	Cleaning kit.

ACCESSORIES FOR OTHER MODELS OF STEREO SERIES

ST-001	Eyepieces (pair) WF5x/22mm.
ST-002	Eyepieces (pair) WF10x/20mm.
ST-003	Eyepieces (pair) WF15x/15mm.
ST-004	Eyepieces (pair) WF20x/13mm.
ST-005	Eyepiece micrometer WF10x/20mm, 10mm/100um.
ST-020	1x objective (for S-10 and S-20 series).
ST-021	3x objective (for S-10 and S-20 series).
ST-022	4x objective (for S-10 and S-20 series).
ST-025	1x objective (for ST-50Led).
ST-026	3.5x objective (for ST-50Led).
ST-040	Darkfield condenser (for all, except ST-45-2L).
ST-041	Sample clip.
ST-032	Eyecups (pair).
ST-011	White/black object-plate, diameter 60 mm (for ST-45-2L).
ST-012	White/black object-plate, diameter 95 mm (for all, except ST-45-2L).
ST-014	Glass stage, diameter 95 mm (for all, except ST-45-2L).
ST-015	Glass stage, diameter 60 mm (for ST-45-2L).
ST-006	Tungsten bulb, 12V/10W.
ST-028	Halogen bulb 12V/10W.
DC-001	Plastic dust cover, small (for all, except ST-50Led).
DC-002	Plastic dust cover, medium (for ST-50Led).
M-114	C-Mount adapter for 1/2" sensor.
M-115	C-Mount adapter for 1/3" sensor.
M-116	C-Mount adapter for 2/3" sensor.
M-173	Photo adapter for APS-C and Full Frame Reflex cameras.
M-113.2	Ring adapter, 30.5mm (for monocular and binocular microscopes).
M-069	Solar battery pack (for ST-30-2Led, ST-30-2LedR, ST-50Led).
M-005	Micrometer slide for software calibration, 1mm/10um, 10mm/100um.
15104	Cleaning kit.

ACCESSORIES FOR LAB SERIES

ST-001.1	Eyepieces (pair) WF5x/22mm.
ST-401	Eyepieces (pair) WF10x/20 mm.
ST-402	Eyepieces (pair) WF15x/15 mm.
ST-403	Eyepieces (pair) WF20x/10 mm.
ST-405	Eyepiece micrometer WF10x/20mm, 10mm/100um.
ST-085	Additional lens 0.5x (w.d. 165mm) (Only for LAB-20 and LAB-30.)
ST-086	Additional lens 1.5x (w.d. 47mm) (Only for LAB-20 and LAB-30.)
ST-415	Eyecups (pair).
DC-002	Plastic dust cover, medium.
ST-092	Protective glass for stereohead.
ST-417	White/black stage-plate, 60mm diameter for LAB Series.
M-114	C-Mount adapter for 1/2" sensor.
M-115	C-Mount adapter for 1/3" sensor.
M-116	C-Mount adapter for 2/3" sensor.
M-173	Photo adapter for APS-C and Full Frame Reflex cameras.
M-113.1	Ring adapter, 30mm (for monocular and binocular microscopes).
ST-418	Focusable C-Mount adapter for 1/3" sensor.
ST-419	Focusable C-Mount adapter for 1/2" sensor.
M-005	Micrometer slide for software calibration, 1mm/10um, 10mm/100um.
15104	Cleaning kit.
VP-LAB	IQ/OQ/PQ Validation Protocols.

M-069 - Solar battery pack

Included battery: rechargeable – Lithium-Poly. Capacity: 2600 mAh.
Output voltage: 5,5Vdc. - Dimensions: 120x73x10 mm.
Autonomy: over 6 hours at medium intensity (X-LED³).
Charging modes: with solar panel (12h),
with external USB power supply (not included)
or from PC USB port (5h).



15104 - Cleaning kit

It cleans glass quickly and effectively,
without leaving residue or odor.
Ideal for precision lens or prism cleaning.



ST-041 - Sample clip



Multi Lab SRL

Str. Baia de Arama, nr.1, Bloc C
Spatiu Comercial Nr. 2,
Sector 2 - Cod.Postal 022204
Bucuresti - Romania

Tel.+4 021 3127144 / 45 :: Fax.+4 021 3127141

www.multilab.ro :: www.labshop.ro :: www.eulab.eu